

less sensitive areas of the body, such as the chest, arms, legs and feet. Use the Precision tip (#150) for tiny and sensitive areas of skin, such as nose.

## HOW TO USE THE NuAcquaDerm™ DEVICE

1. Clean the skin of oil, and make-up. Dry thoroughly.
2. Choose between the diamond heads (fine, coarse or precision) and attach to the handheld unit.
3. Push the ON button/Low suction to begin treatment. To choose High suction push ON button second time.
4. Hold the skin tight with one hand and gently glide the Diamond Head away in the opposite direction. There is no need to press hard on the skin's surface.
5. Keep the Diamond Head flush with the surface of the skin to ensure proper suction. Limit the number of strokes over the same part of the skin to 2-3 passes until your skin becomes used to the treatment.
6. Do not hold the diamond tip in one spot while suction is on.
7. The NuAcquaDerm unit is waterproof, but couldn't be submerged in water or any form of liquid.

## POST TREATMENT CARE

1. Some redness of the skin is normal, and will fade after your treatment.
2. Cleanse your face immediately after treatment with a wet cloth and rehydrating toner to remove any dead skin cells left behind.
3. Once the skin has dried completely, apply a rich moisturizer.
4. Do not use any other method of mechanical or chemical exfoliation for two to five days.
5. Don't forget to use sunscreen every time you go out.
6. Avoid staying in the sun for 7 days after microdermabrasion treatment to prevent UV rays from damaging your skin.

## CLEANING AND MAINTENANCE

Remove the collection cup/container and clean with an antibacterial soap.

Remove diamond head from device. Clean the diamond head with rubbing alcohol (isopropyl).



Dry the diamond head with a cloth and place back on to the device.

Diamond heads are medical grade and recommended every 60 days replacement.

## CHARGING THE UNIT

Prior to first time use, charge Your NuAcquaDerm tool for 8-10 hours – this is for first time use only and to make sure your hand held tool is fully charged. A fully charged NuAcquaDerm gives 45 minutes of use.



## WARRANTY

Your NuAcquaDerm device is warranted to be free from defect in material and workmanship for a period of two (2) years under normal use after its original purchase date.

*If the product should become defective within the warranty period, contact Customer Service at [www.nuacquaderm.com](http://www.nuacquaderm.com)*

This warranty does not cover products damaged by the following:

- Accident, misuse, abuse or alteration.
- Servicing by unauthorized persons.
- Use with unauthorized accessories.
- Any other conditions beyond our control.

For issues or questions, please contact us at:  
888-637-7958  
Monday-Friday, 9:00am - 5:00 pm EST

**NuAcquaDerm hereby disclaims all representations and warranties express or implied, as to the products, including,**

**but not limited to, any implied warranty of merchantability, non-infringement or fitness for a particular purpose. The NuAcquaDerm will not be liable for any damage, loss, cost or expense for breach of warranty.**



# NUACQUADERM

**Hand-held waterproof microdermabrasion beauty device**



## NuAcquaDerm™

Congratulations on your purchase of NuAcquaDerm™ the only filters free at-home waterproof microdermabrasion device of both the face and the body! The NuAcquaDerm is a hand-held, spa-like beauty treatment that sweeps away dead skin cells leaving behind a smoother, healthier-looking and more luminous complexion. Get more beautiful-looking skin – from head to toe!

Please read this Instruction Manual carefully before getting started.

### WHAT IS MICRODERMABRASION?

As we age, the skin cell turnover rate slows down, often leaving the skin looking dull and rough. Microdermabrasion is the process of abrasion and suction. As the diamond tip surface passes over the skin, the device vacuums up the exfoliated particles and skin debris. Through the combination of abrasion and suction, the skin dermis is cleaned and brightened. With regular use, your skin will look thicker, smoother and healthier.

### PRODUCT KEY:

1. Body
2. Collection cup/container
3. Wireless charging base
4. Coarse-surface tip for body #75
5. Coarse-surface tip for face #100
6. Fine-surface tip for face #200
7. Precision tip #150

### NuAcquaDerm™ KIT

- Waterproof, Cordless and Rechargeable NuAcquaDerm tool
- Fine-surface tip for face #200
- Coarse-surface tip for face #100
- Coarse-surface tip for body #75
- Precision tip #150
- AC 100-240 V adapter UL/Universal
- User manual
- 2-years limited warranty



### OPERATING FEATURES

1. Button On (low/high intensity) Off
2. Diamond heads
3. Base o-ring
4. Base plate
5. Light Indicator

### CAUTIONS

If you are under a doctor's care for ANY skin condition, or if you have or suspect that you may have a skin condition, please consult your doctor before using the NuAcquaDerm.

The NuAcquaDerm is recommended for once-a-week use only.

The NuAcquaDerm tool has been designed for professional and personal use. We recommend that you do not share your NuAcquaDerm device.



### WARNINGS

The appliance must only be used on healthy skin.

DO NOT use the NuAcquaDerm if you have the following conditions:

- Skin Cancer
- Herpes outbreak
- Active rosacea (eczema, psoriasis or dermatitis)
- Sunburn, open wounds or sensitive skin
- Recent cosmetic surgery or procedure (consult your doctor)
- Pregnant (consult your doctor)
- DO NOT use any chemicals products such as alpha and beta hydroxy acid, or retinol products prior to or within 48 hours of using the NuAcquaDerm
- DO NOT continue to use the NuAcquaDerm if you feel any pain during the treatment
- DO NOT continue to use the NuAcquaDerm if you experience any adverse effects (rash, swelling, burning, bruising, irritation or discomfort)

### IMPORTANT WARNINGS & CAUTIONS

As with all electrical devices, basic safety precautions should always be followed, especially when children are present. Keep out of the reach of children. Please read all instructions before use.

Important Warnings:

- Use the NuAcquaDerm device only as directed in this instruction manual.
- Only use with the NuAcquaDerm tool attachments.
- Do not use, place or store the NuAcquaDerm charging cradle where it can fall or be pulled into a tub or sink.
- Do not force the AC adaptor into an electrical outlet. This can damage the AC adaptor or the electrical outlet, and may also result in serious injury.
- Do not use the AC adaptor if the cord or the plug is damaged.
- Do not bend or abuse the wire of the power cord and keep away from heated surfaces.
- Do not put the NuAcquaDerm device in your eyes or use on the eyelid.

- Do not shake, drop or open the NuAcquaDerm device. Doing so may void the warranty.
- Always keep the NuAcquaDerm device away from heat over 80° F, humidity over 90% and direct sunlight.
- Do not use the NuAcquaDerm device on children.
- Stop using the NuAcquaDerm device immediately if you think there may be something defective or wrong with it.

### STAINLESS STEEL DIAMOND TIPS

NuAcquaDerm uses professional grade Stainless steel diamond tips as the exfoliating agent of the Microdermabrasion procedure. The NuAcquaDerm kit comes with 4 Diamond Tips of different intensities:

1. Fine-surface tip for Face #200
2. Coarse-surface tip for Face #100
3. Coarse-surface tip for Body #75
4. Precision tip #150

Use the Fine tip (#200) for sensitive areas of skin, such as the face and neck. Use the Coarse tips (#100 and #75) for larger,